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Smart Bracelet User Manual

1. Compatible System of Phone and BLE Version :

Android 4.2 and above


IOS 8.0 and above

BT4.0

2. User Preparation

Please make sure the device is full charged before the first time using.

3. APP " LEFUN HEALTH" Down-load Scan the QR Code

Can download the app " LEFUN HEALTH" from the APPLE store OR Android market

4. Bracelet Connection

First Long press the touch button to turn on the bracelet. Make sure the smart phone bluetooth is on. Second, Open the "LEFUN HEALTH" APP, Press the top left corner, enter into the APP FUNCTION LIST, then press the APP icon--press "Search" choose the right bracelet MAC number will show "Bluetooth Pairing request" press "Pair" .

5. Bracelet Operation

When the bracelet successfully connected with the phone, the bracelet will sync the phone time date. The APP will sync the bracelet sports data, heart rate monitoring data, blood pressure data, and etc. Short press the bracelet touch button to switch the function, long press the touch button to turn ON/OFF.

6. Function Instruction

Home Display: Show Time, Date, Battery and etc;


Pedometer : 24 hours recording on bracelet. Clear bracelet data at 24:00, can check data from the APP;

Sports mode: Statistics of your day Exercise duration, 24:00 Am per day. View history in APP;

Calories/Mileage : 24 hour recording on bracelet. Clear bracelet data at 24:00, can check data from the APP;

Heart rate monitor : Automatically measure the heart rate when the function switch to the Heart rate, or can press the heart rate monitor in APP. Data will be saved in the APP.

Sleep monitoring: Statistics of your day Sleep long clear at 24:00 Am. View history in APP;

Blood pressure : Automatically measure the blood pressure when the function switch to the Heart rate, or can press the heart rate monitor in APP. Data will be saved in the APP.

Smart Alarm: Find the "Smart Clock" on APP, set the time, press the time, can set the repeat mode, press complete to save the setting, the bracelet will vibrate when time arrived.

Sleep monitor : Automatically record your sleep status with analyzing the deep sleep and light sleep hours, data saved in app,

Sedentary Reminder : Open the function in APP and set the time duration, the bracelet will vibrate.

Reminder for drinking : Open the function in APP and set the time duration, the bracelet will vibrate.

Anti-lost : Open the function in APP, the bracelet will vibrate when the phone is taken away.

Message Reminder : Open the function in APP, the bracelet will vibrate when the phone comes a new message (SMS, QQ, Facebook, Twitter, Whatsapp, Wechat.)

Call Reminder : Open the function in APP, the bracelet will vibrate when the phone comes.

Shake the bracelet to take photo : Open the function in APP, then press the "Shake to take a photo", enter into the camera, shake the bracelet to confirm. (Make sure the mobile phone setting: Allow the APP to take photo and save the photo)

Find Bracelet : Press the function in APP, the bracelet will vibrate.

Restore factory Settings: set the time to the previous day in the mobile APP or the data of the bracelet can be cleared by synchronizing one day later;

*The company reserves the right to make changes to the contents of the specification without notice.