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

EN
inskam-306

Important Note

Thank you for purchasing a LCD Digital Microscope. Please take a moment to read through this manual. Its contents will help you to get the best use of the smart microscope.

Before using the LCD Digital Microscope, please ensure that you read and understand the safety precautions described below. Always ensure that the LCD Digital Microscope is operated correctly.

- Keep LCD Digital Microscope away from water and other liquids. Do not use LCD Digital Microscope with wet hands. Never use LCD Digital Microscope in the rain or snow. Moisture creates the danger of fire and electric shock.
- If you notice smoke or a strange odor coming from LCD Digital Microscope, turn the power off immediately.
- LCD Digital Microscope may become slightly warm during use. This does not indicate malfunction.
- If the exterior of LCD Digital Microscope needs cleaning, wipe it with a soft, dry cloth.

Keep away from liquids Stay away from the flame Please use dry cloth for cleaning


EN-1

1 Preparation

(1) Standard accessories

- LCD Digital Microscope
- USB cable
- Cleaning fabric
- Operation guide

(2) Parts Names & Functions

EN-2

11 12 13

- 4.3 inch LCD display
- MicroTF card slot
- HDMI interface
- LED brightness dimmer
- Focusing wheel (Wheel B)
- Microscope tube
- Wheel A
- Object stage
- Capture button
- Reset hole
- DC 5V in
- Charging indicator

(3) Function introduction

- MicroTF card slot: Support MicroTF card (not equipped with this unit) up to 32GB, including FAT, FAT32 two file formats.
- Micro usb: DC 5V power supply and data transmission.
- HDMI high-definition interface: HDMI video transmission.

EN-3


- Power on/off: Short press the power button to enter the file playback mode, long press the switch machine.
- Photograph: Short press to take a photo, long press the video, then press once to stop recording.
- Menu Setup: Short press to enter the menu page, here you can set the Language / LCD Bright / Date and Time / Screen Saver / Power frequency Format / Default Setting / FW Version.
- Zoom in/Right: Photo mode, short press to zoom in, other modes, move up or left.
- Zoom out/Left: Photo mode, short press to zoom out, other modes, move down or right.
- Confirm: Short press to confirm, long press to rotate.
- LED brightness adjustment: Enhance light upwards, weaken light downwards.

EN-4

2 Steps for usage

(1) Standard mode


- Insert the MicroTF card (not equipped with this unit) into the MicroTF card slot; press gently until the card is locked! (Note: MicroTF cards must be inserted to record and store pictures and videos.)
- Long press on/off button, turn on the device.
- Put object to be observed onto the object stage, Rotate Wheel A to bring the object near to the microscope tube.
- Rotate the wheel B left and right to obtain fine focus.
- When a clear focus got, rotate Wheel C clockwise to the object stage.

EN-5

(2) U disk storage mode

Connect the microscope to your computer via a USB cable, turn on the device, select the "mass storage" mode, and view the stored photos and videos on your computer.

EN-6

(3) PC mode

Note: The device button function is no longer active when working in PC mode.

Windows system :

A. Download and install the pieces from the link below www.inskam.com/download/PC.zip

(The browser used is for reference only, you can do it according to your existing browser)

B. Turn on the device and connect the device and computer terminal through the data cable, select PC Camera mode.

C. Open "Smart Camera", click on the device in the settings, as shown below.

(Specific parameter settings can refer to the picture content)

EN-7

Mac system :

A. Turn on the device and connect the device and computer terminal through the data cable, select PC Camera mode.

B. Open the "Photo Booth" software in the Apple App.


C. Click on the Camera and switch the camera as shown below.

EN-8

(4) HDMI mode

Connect the microscope and HD display with HDMI function through HDMI HD cable (not equipped with this unit), then you can observe the images you see in real time on a high-definition display.

EN-8

3 Charge Battery

- Plug the small end of the power adaptor into the DC 5V-in port located on the bottom side of the unit. Then plug the other end of the power adaptor into a standard AC outlet. The battery indicator LED at the end corner of the unit will be in red to show that the battery is charging.
- When the battery is finished charging, the battery indicator LED will be in off. Unplug the unit. The entire charging process takes about 4 hours.
- You can use your microscope for up to 4 hours of continuous use after full charge. When the battery is low, there will be low-power indicator on the left down corner of the panel, red LED light will become weak and image quality becomes poor. Then you have to connect the power adaptor and plug it in to recharge the battery.

EN-9

4 Specifications

Image sensor	2M high quality Image Sensor
Magnification	50x-1000x
Display	4.3 inch TFT panel, 16:9
Focusing range	0-40mm
Focus mode	Manual
Resolution	2048*1536, 1920*1080, 1280*960, 640*480
Data Format	JPG/AVI
Light source	8 LED lights (adjustable brightness)
Frame rate	30fps
Operating system	Windows XP / Vista / Win7/8/10, Mac 10.5 or above
Interface	USB2.0
Lens structure	2G+IR
Aperture	F4.5
Lens angle	16°
Language	English / German / French / Spanish / Italian / Japanese / Simplified Chinese / Traditional Chinese / Portuguese / Russian / OSD language
Power supply	Lithium ion 3.7V / 1000mAh battery Working time: 4 hours; Charging time: 4 hours
Size	300mm x 140mm x 225mm

EN-10

5 Troubleshooting

Problems	Solution
The unit does not turn on	1. Does the battery have electricity? 2. Is the power adapter connected?
Device crash	1. Use a fine pin to insert the back reset hole, reset the device.
The screen is blank	1. Confirm that the lens is be enlarged a correctly position before the shot 2. Whether the focal length has been adjusted to the optimal position
There are smudges or scratches on the screen	1. Clean the screen and lens carefully with soft fabric.
The microscope is connected to the PC and the TV card cannot be read	1. Make sure the HD screen is turned on 2. Whether to successfully switch to the corresponding video input channel on the HD screen 3. Make sure the video cable is securely connected to the microscope video input port
The microscope is connected to the PC, and after downloading and opening the software, there is no image	1. Reconnect the microscope to the PC 2. Confirm whether to enter PC camera mode
The microscope is connected to the PC and the TV card cannot be read	1. Reconnect the microscope and PC 2. Confirm whether to enter the Mass Storage mode
The battery indicator LED is not in red, even though the power adapter is connected and plugged into a wall socket, the battery is not charging	1. Plugged into the Microscope's DC-in port 2. Ensure that the outlet has power and is not malfunctioning. If the outlet is connected to a wall switch, make sure that the switch is on

EN-11