
H1

Health Smart Watch

Quick start manual


H1

- EN** *Quick start manual*
- DE** *Einführungsanleitung*
- FR** *Manuel de démarrage rapide*
- JP** *クイックスタートガイド*
- ES** *Manual de Inicio Rápido*
- IT** *Manual di istruzioni*

H1


EN

Quick start manual

1. Package Contents

EN

Watch	Charging Cable	Quick Start Manual
		

2. How to Operate


2.1 Powering on Your Watch

Long press the power key to turn on your watch. Please make sure to fully charge the watch first if it Doesn' t turn on.

2.2 Charging Your Watch

EN

Properly attach the magnetic charging base to the back of the watch, then plug the USB end into a USB port on the computer, charging dock or power bank for charging.

Input current: $<0.3A$

Input voltage: DC5V


Charging time: about 2 hours

Note:

It is recommended to use the universal charger with certification mark on the market.

2.3 Wearing Your Watch

1. For optimized tracking of the measured values, we recommend wearing the watch with one finger's width below your wrist bone.
2. Please make sure the watch is fairly snug against your skin and does not slide up or down your wrist during exercise.

2.4 Navigating Your Watch

EN


1.Key Instruction

	Short Press	Long Press
Power Key	1.Return to the previous menu 2.Pause / Continue the exercise	1.Power on or off 2.To end the exercise from the activity pausing interface

2.Gestures Instruction

Tap the screen	Confirm to use this feature/Enter the sub-interface
Swipe up/down	1.Slide down the home screen to view the workout and heart rate data today 2.Switch screen

Dial Frame

Indicator Lamp

3.Icon Instruction

EN

Bluetooth connected	Low Power	Charging	Fully Charged	End of Workout	Continue Workout
					

2.5 APP Downloading

Download the “Zeroner Health Pro” app from the App Store or Google Play Store, or scan the QR code below to download.

2.6 Pairing Your Watch with Your Phone

After downloading the APP, please follow the following operating instructions to complete the connection between the watch and the mobile phone:

- 1.To successfully pair your watch and mobile phone, you need to turn on Bluetooth of your mobile phone;
- 2.Open APP, register an account, and truthfully fill in your personal information (height, weight, date of birth);
- 3.Click "Device" > "Watch+" > "H1-****" (Watch Model) > Connected successfully.

Notes:


- 1.The app needs permission to access your location if you are using an Android phone, so you need to **ALLOW LOCATION ACCESS** in System Settings in order to connect successfully.
- 2.The first time you connect to the app, the date & time on your phone will be synchronized to the watch, and previous data for steps, distance & calories on the watch will be cleared.

2.7 Replacing the Strap

EN

Please choose the strap with a width of 20mm if you want to replace it.

1. Remove the strap from the watch by sliding the snap lock on the strap.
2. Align the new strap with the watch and buckle the strap in.

3.Main Features

EN

3.1 Heart Rate Monitoring

- 1.The watch can monitor 24-hour real-time heart rate, resting heart rate and heart rate zones during exercise, you can choose to enable or disable this function on the APP.
- 2.Real-time heart rate measurement: Wear the watch correctly, and down slide the watch from the home screen to the interface of heart rate measurement. Please wait for a few seconds with the watch screen upward, and wait for additional a few seconds after obtaining the value to make the result more accurate.
- 3.Multiple heart rate zones display during exercise: all detailed data can be viewed after connecting and syncing to the app.

Note:

The light signal transmission might be blocked if your skin is too dark or with excess hair, or improper wearing also may lead to measurement failure.

3.2 ECG

EN

1. Wear the watch correctly, slide upward to the “ECG ” feature interface from the main screen, and click the icon to start measurement; lightly press your thumb on the dial frame on one side of the watch, fix your index finger on the opposite side, and wait for the measurement progress box to finish. Please be careful not to touch your right hand with your left hand.
2. After the measurement, if you need to play back the ECG, please view it after synchronizing your data on the APP. The measurement results are for reference only, and can not be used as a basis for clinical diagnosis.

Notes:

1. **Please relax your fingers when measuring, and do not exert too much force, otherwise the measurement result will be affected.**
2. **If your skin is too dry or too hairy when measuring, the electrical signal transmission might be affected, which may lead to measurement failure.**

3.3 Vascular Health (APG)

EN

1. Wear the watch correctly, slide upward to the “APG ” feature interface from the main screen, and click the icon to start measurement. Please press your thumb lightly on the Indicator Lamp area, keep sitting still and wait for the measurement progress box to finish.
2. For more information, please view if after synchronizing your data with the watch on the APP. The measurement results are for reference only, and can not be used as a basis for clinical diagnosis.

3.4 Blood Pressure and Blood Oxygen (BP&SpO2)

1. Upon the first measurement of blood pressure or after the factory resetting of the watch, please use the home sphygmomanometer for blood pressure calibration on the APP.

Calibration Steps: APP> ‘Device’ > ‘Blood Pressure Calibration’ , then follow the on-screen instructions. According to the statistical data, the more times the watch is calibrated, the closer the results will be to the value of the sphygmoma-

nometer. Please keep a home sphygmomanometer on hand and develop a good habit of regular calibration.

2. Please rest for 10 minutes before the blood pressure measurement, then put on the watch properly, put your arm flat on the support, keep sitting still and relax your whole body naturally.

3. Slide upward to the "BP&SpO2" from the main screen, and click the icon to start measurement. Please press your thumb lightly on the Indicator Lamp area, and wait for the measurement progress box to finish.

4. After the measurement, the results can be read in real time on the watch. For detailed information, please view it after synchronizing your data with the watch on the APP. The results are for reference only, and can not be used as a basis for clinical diagnosis.

Notes:

EN

1. There should be no movement during the blood pressure measurement, otherwise the false wave due to muscle movement may cause the measurement failure.
2. It is recommended that the following populations should concern more about the blood oxygen saturation SpO₂: the high-intensity brain workers, the aged population, and people with hypoxia in the high mountain/highland environments.

3.5 Heart Health and Scientific Sleep

1. Wear the watch correctly, slide upward to the “Heart Health” feature interface from the main screen. Turn this function on or off as needed.
2. When this function is turned on, the watch will automatically monitor the heart health at night. If you wear the watch for more than 6 hours, you can view a scatter chart of your heart health after synchronizing your data with the watch on the APP. The measurement results are for reference only, and can not be used as a basis for clinical diagnosis.

3. When this function is turned on, the watch will enter the sleep monitoring state from 8:00 pm to 9:00 am the next day. You can check REM sleep and other sleep data after synchronizing your data with the watch on the APP, so as to know about your sleep status.

Notes:

- 1. When this function is turned on, the heart rate lamp will always be on and the device battery life will be shortened to about two days.**
- 2. The sleep monitor will stop recording within 5-10 minutes after you wake up to move.**
- 3. The watch will not record daytime sleep data.**

3.6 Temperature

1. Wear the watch correctly, slide upward to the "Body Temp" feature interface from the main screen, and click the icon to start measurement. Please turn the watch screen upward while measuring and wait for the measurement progress box to

finish. Wait for a few seconds after getting the value, so that the result will be more accurate.

2. Once the measurement is complete, the watch will display the skin temperature and the estimated body temperature in real time. For more information, please view it after synchronizing your data with the watch on the APP.

Note:

When you just start to wear the watch or after a drastic change in the ambient temperature, please wait until the contact temperature between the watch and the skin returns to the normal state before the temperature measurement.

3.7 Body Fat Rate

1. Please make sure the personal information (gender, height, and weight) on the watch is correct before measurement, otherwise the incorrect values will affect the measurement accuracy.

2. Wear the watch correctly, slide upward to the "Body Fat" feature interface from the main screen, and click the icon to start measurement. Lightly press your thumb on the dial frame in measurement, fix your index finger on the opposite side, and wait for the measurement progress box to finish. Please be careful not to touch your right hand with your left hand.
3. After the measurement is completed, the watch will display the results in real time. For more information, please view it after synchronizing your data with the watch on the APP. To improve accuracy, take measurements in the same position and at the same time of the day.

Notes:

1. Try to avoid measurement after strenuous exercise, sauna or shower, excessive alcohol and water intake, or meals (1-2 hours).
2. The following people may have a large difference between the measurement results and the actual body fat rate due to large differences in the average water content or the bone mineral density, etc.: children in the formative years, the elderly, people with fever, people with osteoporosis, pregnant

women, bodybuilders, sports practitioners, etc.

3.Low skin and watch temperature or too dry skin may cause measurement failure.

3.8 Stress Value

1. Wear the watch correctly, slide upward to the “Stress” feature interface from the main screen, and click the icon to start measurement. Please turn the watch screen upward while measuring, keep sitting still and wait for the measurement progress box to finish.

2. After the measurement is completed, the watch will display the results in real time. The higher the value, the greater the pressure. For more information, please view it after synchronizing your data with the watch on the APP.

Note:

It is recommended to measure at the same time period of the day. It is best

measured in the quiet state after getting up every day.

3.9 Sports

1. The watch supports two exercise modes: walking and free training. You may choose as needed, set your goals and get real-time statistics.
2. During the exercise, short press the power button to pause the exercise, and long press the power button to exit the exercise mode after the pause.
3. After the exercise is over, the watch will display the exercise data in real time. For more information, please view it after synchronizing your data with the watch on the APP.

4.1 Water Resistance Instruction

The water resistance rating: IP65

The water resistance performance of the device is not permanently valid, it may diminish as time passes. The device can be used during washing hands, rain, but not support hot water shower, swimming, diving, surfing, etc. It does not have waterproof effect on corrosive liquids such as acidic and alkaline solutions, and chemical reagents. If you inadvertently encounter corrosive liquid, please clean with clear water and wipe it dry. The damage caused by abuse or improper use is not covered by the warranty.

The following conditions may affect the water resistance performance, and should be avoided when using:

1. The watch falls, bumps, or suffers from other impacts.
2. The watch is exposed to detergent, lotion, oil etc.
3. Hot and humid scenes such as hot baths and saunas.

4.2 Specification Parameter

EN

Physical size	45*34.5*10.5mm	Adjustable strap	160mm-210mm
Display size	1.3 inch TFT color square display	Working temperature	-10-45 °C
Weight	About 20g (with strap)	Screen resolution	240*240 pixels
Battery capacity	170mAh Li-Polymer battery		
Battery Life	5-7days (Receive an average of 50 messages & 5 calls per day;raise your wrist to wake screen up 50 times;exercise for half an hourper day;turn on24-hour automatic heart rate monitoring) .		

4.3 Firmware Upgrade

1.Firmware Upgrade

When there is an new firmware version, there will be a notification in the app. Go to the "Device" interface of the app and select firmware upgrade.

Notes:

- 1. Make sure that the battery level is more than 50% before upgrading.**
- 2. During the upgrade process, keep the screen of your phone bright, and only when the upgrade is completed can you exit the interface, otherwise, the upgrade will fail.**

2. Upgrade Failed

Wait for the watch to restart automatically if upgrade fails. Then reconnect your watch to the app for upgrading again.

Note: allow the following conditions to ensure a successful upgrade.

- 1) The Bluetooth of the mobile phone and that of the watch are connected;**
- 2) Smooth access to the network;**
- 3) The APP has been granted the location permission;**
- 4) The APP has been granted the storage permission.**

5. Device maintenance

EN

5.1 Device Care

1. Do not use any sharp objects to clean the device.
2. Avoid using solvents, chemical cleaners or insect repellents that might damage plastic components of the device.
3. Thoroughly rinse the device with clean water after exposure to chlorine, salt water, cosmetics, alcohol or other harsh chemicals to avoid damage to the device.
4. Avoid pressing the key on the device while it is under water.
5. Avoid extreme shock and harsh treatment of the device, as it may reduce the product life.
6. Do not exposure the device to extreme temperatures that are either too cold or too hot for a long time, which may cause permanent damage.
7. After each training session, please rinse the watch with clean water.

5.2 Cleaning the Device

1. Gently wipe the device by using a cotton cloth with a mild neutral detergent.
2. Please allow the device to dry completely after cleaning.

Note:

Even the faintest prolonged exposure to sweat or moisture may cause corrosion of the charging terminal when charging the device, which will also hinder the transmission of data and affect the charging.

5.3 Keep your watch and wrist dry

1. Be sure that your wrist is dry before you put on your watch.
 2. After washing or sweating, dry it thoroughly before putting it back on your wrist.
- Don't wear your watch too tight.


6.Important Safety Tips


- 1.Always consult with your physician before beginning or modifying any exercise program.
- 2.The device, accessories, heart rate monitor and related data are intended to be used only for exercise monitoring not medical purposes.
- 3.The heart rate readings are for reference only, and no liability is accepted for the consequences of any wrong interpretation.
- 4.Do not expose the watch to a heat source or in a high-temperature location, for example, in an unattended vehicle in the sun. To avoid the possibility of damage, remove the device from the vehicle or store it out of direct sunlight.
- 5.If you want to store the watch for an extended time period, please store it within the temperature ranges specified in this manual.

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